



## I.V. Flow regulator - Dosy 2007

Version with 4 indented grips in the lower part and female tube connector

**Code 8008**

### INTENDED USE

Regulation of fluid delivery in I.V. administrations, for gravity use only.

### MAIN FEATURES

This device ensures constant administration over time. It allows extreme ease and precision in setting the flows.



### FEATURES

#### Body made with high-heat resistant ABS

This type of material allows the maintenance of mechanical characteristics over time.

#### Gasket made in clear type silicone LSR

Silicone is the best material for the gasket, because it combines exceptional mechanical characteristics, a perfect biocompatibility and the absence of interaction with drugs.

#### Microchannel patented

The distribution channel is designed with a specific shape that makes delivery much more accurate, especially at low rates.

#### Easy and smooth regulation

The silicone gasket is characterized by a low coefficient of friction which allows a smooth rotation of the piece.

#### Low standard deviation, high precision



The accuracy of the processing and the solutions adopted guarantee a high uniformity of delivery.

#### Ink type

Scales and logo are printed with UV cured ink to achieve a good stability towards both ethyl alcohol and mechanical stress.

#### Graduated scale

The flow regulator, sometimes referred to as dial-a-flow, is available with single scale, or other scales. Possibility to print the customer's logo. Several colours are available.



## SHELF LIFE

The 5-year shelf life has been verified by an external laboratory with a study of natural ageing on samples kept in correct condition in warehouse for five years.

The product was found to be compliant after ageing.

## DELIVERY RANGE

5-300 ml/h

## SCALE ACCURACY

Samples are tested at flow rates of 20, 50, 60, 80, 125 and 200 ml/h and the delivered quantities are recorded after 60 minutes, with a  $\Delta h$  of 100 cm.

The accepted deviation from the nominal value at 20 ml/h is  $\pm 20\%$ .

The accepted deviation from the nominal value at 50, 60, 80, 125 and 200 ml/h is  $\pm 15\%$ .

## CONSTANCY OF DELIVERY

Samples are tested at a flow rate of 50 ml/h and the delivered quantities are recorded each hour over a period of 6 hours.

The accepted deviation from the mean value is  $\pm 10\%$ .

## HEAT STABILITY - NO LEAKAGE

Dosy 2007 is made of heat resistant materials, which guarantee stability after thermal treatment and ageing.

There is no evidence of leakage for samples that were treated at 65 °C for 7 days. Sealing tests according to ISO 8536-13:2016 at 0.5 bar.

Tested in OFF and in OPEN position, with air and water.

## TUBING PORT

Female connector suitable for assembling with tubing.

Inner  $\varnothing$  (ID) 2.74 mm (0.108")

Outer  $\varnothing$  (OD) 3.86 mm (0.152")

## STANDARD

Dosy 2007 complies with the following standards.

### ISO 8536-13:2016

Infusion equipment for medical use.  
Part 13: graduated flow regulators for single use with fluid contact.

### ISO 8536-4:2019

Infusion equipment for medical use  
Part 4: infusion sets for single use, gravity feed.

## BIOLOGICAL EVALUATION

Tested according to ISO 10993-1

## STERILIZATION

EtO - temperature up to 60°C

Gamma - 25 kGy



Phoenix R&D s.r.l.  
Viale della Cooperazione, 9 • 45100 Rovigo (RO) - Italy  
tel. +39 0425 471507 • info@phoenix-rd.com  
www.phoenix-rd.com

Regional Office for UK  
info.uk@phoenix-rd.com  
mob. +44 7941565386



QMS Cert nr. 501007251